



User's Manual

PTEN Activity ELISA Kit



DEIA-XYZ23



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The PTEN Activity ELISA is designed to quantify the phosphatase activity of PTEN by detection of the product, PI(4,5)P₂, in a competitive ELISA format, eliminating the need for radioactivity, organic solvents, and thin layer chromatography.

General Description

PTEN (Phosphatase and Tensin Homolog deleted on Chromosome 10) is a 3' phosphoinositide phosphatase that converts PI(3,4,5)P₃ to PI(4,5)P₂, thus opposing PKB/Akt activation by PI 3-K. PTEN is involved in neuronal stem cell proliferation and self-renewal, cardiac myocyte hypertrophy and contractility, and a wide range of developmental processes. PTEN, however, is best known for its role as a tumor suppressor. Loss of PTEN activity results in accumulation of PI(3,4,5)P₃, abnormal activation of PKB/Akt, unregulated cell growth, suppression of apoptosis, and increased tumorigenesis in a number of human tissues. It has also been proposed that PTEN is a candidate for targeted chemotherapy because certain anti-cancer agents preferentially destroy tumors with PTEN mutations. In addition to this direct role in cancer, PTEN has recently been shown to regulate cancer-associated pathways including VEGF-mediated angiogenesis and others.

The PTEN Activity ELISA directly detects the phosphoinositide product compared to other assays, which detect free phosphate. This eliminates many possible sources of error due to the fact that inorganic phosphate is the product of many phosphatase enzyme activities, and is found in common buffers and cleaning products.

Principles of Testing

The assay is a competitive ELISA in which the signal is inversely proportional to the amount of PI(4,5)P₂ produced. After the PTEN reactions are complete, reaction products are added to the PI(4,5)P₂-coated microplate and a PI(4,5)P₂ detector protein is then added for competitive binding. A peroxidase-linked secondary detector and colorimetric detection is used to detect PI(4,5)P₂ detector binding to the plate. The colorimetric signal is inversely proportional to the amount of PI(4,5)P₂ produced by PTEN.

Reagents And Materials Provided

1. PI(4,5)P₂ Coated Strip-well Detection Plate: 1 plate
2. PI(3,4,5)P₃ Substrate, diC16: 30 nmol
3. PI(4,5)P₂ Standard, diC16: 12 nmol
4. PI(4,5)P₂ Detector: 2 x 2.5 µg
5. 5x PTEN Reaction Buffer: 4 mL
6. Protein Stabilizer: 600 µL
7. Secondary Detector: 300 µL
8. PBS Tablet: 1 tablet
9. 10x PBS-T Buffer: 30 mL

10. TMB Solution: 12 mL
11. 1 N H₂SO₄ Stop Solution: 10 mL
12. DTT: 3 x 50 µmol
13. Yellow 96-well polypropylene U-bottom plate: 1 plate
14. Clear acetate sheet, 1 side adhesive: 3 seals

Materials Required But Not Supplied

1. Source of PTEN Enzyme.
2. Incubated plate shaker or 37°C incubator.
3. Microplate Reader with capability to read absorbance at 450 nm.

Storage

The kit comes in two parts with different storage requirements. Upon receipt store Kit Part 1 at 4°C and Kit Part 2 at -20°C. Store prepared reagents as indicated in the protocol.

Assay Procedure

Please read through entire protocol carefully along with the "Assay Notes" before beginning the assay.

Reagent Preparation

1. PTEN Reaction Buffer:

Prepare fresh PTEN Reaction Buffer and DTT for use on the day of the assay. Dilute the 5x PTEN Reaction Buffer 5-fold in dH₂O and supplement with 10 mM DTT. Each vial of K-DTT1 contains 50 µmol DTT. Add 50 µL dH₂O for a 1M stock. Once reconstituted, the DTT should be used immediately and any remaining should be discarded.

5 mL of PTEN Reaction Buffer = 1 mL of 5x PTEN Reaction Buffer + 50 µL 1M DTT + 3,950 µL dH₂O.

2. PBS Buffer:

Prepare the PBS Buffer by dissolving the provided PBS tablet in 200 mL dH₂O.

3. PBS-T Buffer:

Prepare the 1x PBS-T Buffer by diluting 30 mL of the 10x PBS-T Buffer with 270 mL dH₂O.

4. PI(3,4,5)P₃ Substrate:

Equilibrate vial of PI(3,4,5)P₃ Substrate to room temperature. Prepare a 100 µM PI(3,4,5)P₃ Substrate stock solution by adding 300 µL dH₂O to the vial of PI(3,4,5)P₃ Substrate. Vortex for at least 60 seconds to resuspend the lipid. Spin down and place vial at room temperature. Prior to use, dilute the required amount of PI(3,4,5)P₃ Substrate in PTEN Reaction Buffer for a 16 µM (2x conc.) working solution. 1 mL of 16 µM PI(3,4,5)P₃ Substrate can set up 30 PTEN reactions. The unused portion of 100 µM PI(3,4,5)P₃ Substrate stock can be stored at -20°C for up to 3 months.

- 1 mL of 16 µM PI(3,4,5)P₃ Substrate = 160 µL of 100 µM PI(3,4,5)P₃ Substrate stock + 840 µL of PTEN Reaction Buffer.

5. PI(4,5)P2 Standard:

Equilibrate vial of PI(4,5)P2 Standard to room temperature. Prepare a 40 μ M PI(4,5)P2 Standard stock solution by adding 300 μ L dH₂O to the vial of PI(4,5)P2 Standard. Vortex for at least 60 seconds to resuspend the lipid. Spin down and place vial at room temperature. The unused portion of 40 μ M PI(4,5)P2 Standard stock can be stored at -20°C for up to 3 months.

PTEN Reaction and Incubation

1. Isolate or prepare PTEN according to usual protocols. See attached support protocol for immunoprecipitation of PTEN from cells. Prior to use, dilute the required amount of PTEN to a 2x concentration in the PTEN Reaction Buffer.
2. PTEN reactions can be set up in micro centrifuge tubes or in the yellow Incubation Plate (using recombinant PTEN enzyme only). Please read assay notes at the end of the protocol first.
 - a. For each 60 μ L PTEN reaction (for duplicate assay points): combine 30 μ L of the 16 μ M PI(3,4,5)P3 Substrate (480 pmol) and 30 μ L of the PTEN (2x conc.). You may also include an Enzyme Only control by replacing 30 μ L of Substrate with 30 μ L of PTEN Reaction Buffer.
 - b. Seal the PTEN reactions and let them proceed without vigorous shaking for a certain time, usually 1-4 hours at 37°C or RT Note: If using PTEN bound to beads, light agitation or slow rotation will be necessary to keep beads in suspension.
3. Stop each 60 μ L reaction.
 - a. If PTEN enzyme is bound to beads, centrifuge to separate the beads and transfer reaction supernatant to clean tubes.
 - b. For recombinant PTEN, heat the reaction for 3 minutes at 95°C to stop the reaction.
4. Add an additional 62 μ L of PTEN Reaction Buffer to each stopped reaction for a total of 122 μ L.

Note: Reaction products can be stored at -20°C for up to a week. ELISA Detection can be run on another day.

We suggest that standards and controls be run in duplicate or triplicate. The Incubation Plate layout shown on the next page gives an example plate layout with triplicate standards and controls.

5. Prepare PI(4,5)P2 Standards and Controls:

- a. From the 40 μ M PI(4,5)P2 Standard stock prepared earlier, prepare a 4 μ M working solution by adding 40 μ L of the 40 μ M stock solution to 360 μ L of the PTEN Reaction Buffer.
- b. Make four, 2-fold serial dilutions from the 4 μ M PI(4,5)P2 stock with PTEN Reaction Buffer.

e.g. Each 400 μ L dilution = 200 μ L previous dilution + 200 μ L PTEN Reaction Buffer.
- c. Prepare a 4 μ M working solution of PI(3,4,5)P3 Substrate for the "No Enzyme" control by diluting 8 μ L of the 100 μ M PI(3,4,5)P3 Substrate stock solution with 192 μ L PTEN Reaction Buffer.
- d. Pipet 60 μ L/well of each PI(4,5)P2 Standard solution to rows A through E of the Incubation Plate.
- e. Pipet 60 μ L/well of the 4 μ M PI(3,4,5)P3 Substrate solution to the No Enzyme control wells in row F of the Incubation Plate.
- f. Pipet 60 μ L/well of PTEN Reaction Buffer to the No Lipid control wells in row G of the Incubation Plate.
- g. Pipet 120 μ L/well of PTEN Reaction Buffer to the Blank control wells in row H of the Incubation Plate.

	Standard	Reactions		
A. 200 pmol	○	○		
B. 100 pmol	○	○		
C. 50 pmol	○	○		
D. 25 pmol	○	○		
E. 12.5 pmol	○	○		
F. No Enzyme	○	○		
G. No Lipid	○	○		
H. Blank	○	○		

PIP ₂ Standards and Controls	PIP ₂ per 50 μ L solution
4 μ M	200 pmol
2 μ M	100 pmol
1 μ M	50 pmol
0.5 μ M	25 pmol
0.25 μ M	12.5 pmol
No Enzyme control	0 pmol (with Substrate)
No Lipid control	0 pmol

6. Transfer 60 μ L/well of each stopped PTEN reaction into 2 wells of the Incubation Plate for duplicate data points.

7. Dilute the PI(4,5)P₂ Detector to 0.25 μ g/mL in Detection Buffer:

a. Prepare 10 mL of Detection Buffer by adding 200 μ L Protein Stabilizer to 9.8 mL of PBS Buffer.

b. Pipette 1 mL of Detection Buffer into a vial of PI(4,5)P₂ Detector. Mix gently to reconstitute the Detector. Spin down and transfer the 1 mL of reconstituted PI(4,5)P₂ Detector into the 9 mL of Detection Buffer.

Note: Once the PI(4,5)P₂ Detector has been reconstituted in Detection Buffer, it is only good for 1 day. It cannot be stored for later use.

8. Incubate with PI(4,5)P₂ Detector.

a. Add 60 μ L/well of the 0.25 μ g/mL PI(4,5)P₂ Detector to all Control, Standard and PTEN reaction wells except the Blank controls in row H.

b. Seal the Incubation Plate and incubate for 60 minutes at room temperature with gentle agitation on a plate shaker.

Detection

1. Following the incubation, transfer the reacted mixtures to the Detection Plate (clear flat-bottom strip plate). Transfer 100 μ L from each well to the corresponding well in the Detection Plate. (This can easily be accomplished with a multi-channel pipettor.) Seal the plate and incubate for 60 minutes at 37°C with agitation on a plate shaker.

Note: Never let the detection plate dry out after the assay has started. Always get the next solution ready before discarding the current one from wells in use.

2. During the incubation prepare the Secondary Detector. Briefly centrifuge the vial of Secondary Detector). Dilute the Secondary Detector 1:60 with PBS-T. Dilute ONLY the amount you will use for the current assay and store the remainder of the Secondary Detector at 4°C for future use.

3. After incubation, discard solution from the Detection Plate and wash the wells 3 times with 200 μ L/well of PBS-T.

4. Add 100 μ L of diluted Secondary Detector to each well of the Detection Plate. Seal the plate and incubate for another 30 minutes on a plate shaker at room temperature.

5. Discard the Secondary Detector from the Detection Plate and wash the wells 3 times with 200 μ L/well of PBS-T.

6. Immediately add 100 μ L of TMB solution to each well. Allow color to develop for 15-30 minutes in the dark. Watch for blue color development and DO NOT overdevelop. Stop color development by adding 50 μ L of 1 N H₂SO₄ Stop Solution to each well when the color has turned dark ocean blue in the No Lipid control wells but is still clear to very faint blue in the 200 pmol PI(4,5)P₂ Standard wells. Blue color will change to yellow color upon addition of Stop Solution. Eliminate any big air bubbles present in wells before reading the plate.


Caution: Use caution when dealing with corrosive 1 N H₂SO₄ Stop Solution.

7. Read absorbance at 450 nm on a plate reader.

Assay Notes:

1. Optimization of the PTEN reactions by enzyme amount, reaction time and temperature may be required.
2. The provided 5x PTEN Reaction Buffer must be used for the PTEN reactions.
3. Vigorous shaking will kill PTEN enzyme activity. When using PTEN bound to beads, some agitation will be necessary to keep beads in suspension. When using recombinant PTEN there is no need to shake the reaction.
4. PTEN reactions can be carried out at either 37°C or at room temperature. When performing reactions at room temperature, the incubation time should be increased.
5. The amount of enzyme to use per PTEN reaction will vary according to your individual experiment. Whether you are using purified PTEN or enzyme immunoprecipitated from cell lysate, you will need to try reactions using different amounts of enzyme to determine the optimum condition. When using purified recombinant PTEN from CD, enzyme concentration of 0.5-2 ng/μL is suggested as a starting point. In testing, we found that enzyme immunoprecipitated from cell lysate containing 1-5 mg cellular protein is usually sufficient for each reaction.
6. Beads need to be removed from enzyme reactions before reactions are stored or detected.
7. Step 1 of the "Detection" protocol requires an incubation at 37°C with shaking. If an incubated plate shaker is not available, then the incubation can be done at 37°C without shaking for a period of 2 hours.
8. The detection plate is composed of 12 of 8-well strips. Unused strip wells should be removed from the plate frame and stored in a clean sealable plastic bag at 4°C. Save the plate frame after assay for future use of the remaining strip wells.
9. Never let the detection plate dry out after the ELISA assay has started. Always prepare the next solution needed before discarding the current one from wells in use.

Interpretation Of Results



PTEN activity can be estimated by comparing the absorbance values from the wells containing enzyme reaction products to the values in the standard curve. Plot the absorbance values obtained vs. log of PI(4,5)P₂ in pmol per standard to generate a standard curve using sigmoidal dose-response (variable slope) correlation. Determine the PI(4,5)P₂ level in pmol by interpolation from absorbance values obtained from the enzyme reactions. PTEN activity in your samples can be estimated by the percentage conversion from initial 200 pmol of PI(3,4,5)P₃ per assay point.

References

- (1) T. Maehama and J.E. Dixon (1998) The tumor suppressor, PTEN/MMAC1, dephosphorylates the lipid second messenger, phosphatidylinositol 3,4,5-trisphosphate. *J Biol Chem*, 273, 13375-8.
- (2) T. Maehama and J.E. Dixon (1999) PTEN: a tumour suppressor that functions as a phospholipid phosphatase. *Trends Cell Biol*, 9, 125-8.
- (3) M. Tamura, J. Gu, E.H. Danen, T. Takino, S. Miyamoto and K.M. Yamada (1999) PTEN interactions with focal adhesion kinase and suppression of the extracellular matrix-dependent phosphatidylinositol 3-kinase/Akt cell survival pathway. *J Biol Chem*, 274, 20693- 703.
- (4) L.C. Cantley and B.G. Neel (1999) New insights into tumor suppression: PTEN suppresses tumor formation by restraining the phosphoinositide 3-kinase/AKT pathway. *Proc Natl Acad Sci U S A*, 96, 4240-5.
- (5) J.M. Penninger and J. Woodgett (2001) Stem cells. PTEN--coupling tumor suppression to stem cells? *Science*, 294, 2116-8.
- (6) M. Groszer, R. Erickson, D.D. Scripture-Adams, R. Lesche, A. Trumpp, J.A. Zack, H.I. Kornblum, X. Liu and H. Wu (2001) Negative regulation of neural stem/progenitor cell proliferation by the Pten tumor suppressor gene in vivo. *Science*, 294, 2186-9.
- (7) G. Schwartzbauer and J. Robbins (2001) The tumor suppressor gene PTEN can regulate cardiac hypertrophy and survival. *J Biol Chem*, 276, 35786-93.
- (8) M.A. Crackower, G.Y. Oudit, I. Kozieradzki, R. Sarao, H. Sun, T. Sasaki, E. Hirsch, A. Suzuki, T. Shioi, J. Irie-Sasaki, R. Sah, H.Y. Cheng, V.O. Rybin, G. Lembo, L. Fratta, A.J. Oliveira-dos-Santos, J.L. Benovic, C.R. Kahn, S. Izumo, S.F. Steinberg, M.P. Wymann, P.H. Backx and J.M. Penninger (2002) Regulation of myocardial contractility and cell size by distinct PI3K-PTEN signaling pathways. *Cell*, 110, 737-49.
- (9) B. Stiles, M. Groszer, S. Wang, J. Jiao and H. Wu (2004) PTENless means more. *Dev Biol*, 273, 175-84.
- (10) N.R. Leslie and C.P. Downes (2004) PTEN function: how normal cells control it and tumour cells lose it. *Biochem J*, 382, 1-11.
- (11) Furnari, F. B., Huang, H. J., and Cavenee, W. K., The phosphoinositol phosphatase activity of PTEN mediates a serum-sensitive G1 growth arrest in glioma cells, *Cancer Res*, 58, 5002 (1998).
- (12) D. Haas-Kogan, N. Shalev, M. Wong, G. Mills, G. Yount and D. Stokoe (1998) Protein kinase B (PKB/Akt) activity is elevated in glioblastoma cells due to mutation of the tumor suppressor PTEN/MMAC. *Curr Biol*, 8, 1195-8.
- (13) D.H. Teng, R. Hu, H. Lin, T. Davis, D. Iliev, C. Frye, B. Swedlund, K.L. Hansen, V.L. Vinson, K.L. Gumpfer, L. Ellis, A. El-Naggar, M. Frazier, S. Jasser, L.A. Langford, J. Lee, G.B. Mills, M.A. Pershouse, R.E. Pollack, C. Tornos, P. Troncoso, W.K. Yung, G. Fujii, A. Berson, P.A. Steck and et al. (1997) MMAC1/PTEN mutations in primary tumor specimens and tumor cell lines. *Cancer Res*, 57, 5221-5.
- (14) G.B. Mills, Y. Lu and E.C. Kohn (2001) Linking molecular therapeutics to molecular diagnostics: inhibition of the FRAP/RAFT/TOR component of the PI3K pathway preferentially blocks PTEN mutant cells in vitro and in vivo. *Proc Natl Acad Sci U S A*, 98, 10031-3.
- (15) J. Huang and C.D. Kontos (2002) PTEN modulates vascular endothelial growth factor-mediated signaling and angiogenic effects. *J Biol Chem*, 277, 10760-6.
- (16) Evangelia A Papakonstanti, Anne J Ridley, Bart Vanhaesebroeck (2007) The p110d isoform of PI 3-kinase negatively controls RhoA and PTEN, *The EMBO Journal*, 26, 3050–3061



(17) Jürgen Knobloch, Ingo Schmitz,2 Katrin Götz, Klaus Schulze-Osthoff, and Ulrich Rüther (2008) Thalidomide Induces Limb Anomalies by PTEN Stabilization, Akt Suppression, and Stimulation of Caspase-Dependent Cell Death, Molecular and Cellular, 529- 538



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com