




# Florfenicol and Thiamphenicol ELISA KIT



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

## **Creative Diagnostics**

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221

# PRODUCT INFORMATION

#### **Intended Use**

This kit can be used in quantitative and qualitative analysis of florfenicol and thiamphenicol residue in animal tissue (chicken, duck, pork, chicken liver, pig liver), fish and shrimp, honey, egg, milk (raw milk, UHT pure milk) and feed.

# **General Description**

Florfenicol is a derivative of thiamphenicol, which is broadly applied in animal industry for controlling and treating diseases for it has a relative low cost and strong inhibition against a number of gram bacterium and mycoplasma. The residue of this drug will lead to potential aplastic anemia, so the MRL of which is now restricted to 0.1ppm in food by CAC and many other countries.

This kit is a new product for drug residue detection based on ELISA technology, which can considerably minimize operation errors and work intensity compared with instrumental analysis.

# **Principles of Testing**

This kit is based on Competitive ELISA technology. The microtiter wells are coated with coupling antigen. Florfenicol and thiamphenicol residue in the sample competes with the antigen coated on the microtiter plate for the antibody. After the addition of enzyme labeled anti-antibody, TMB substrate is used to show the color. Absorbance of the sample is negatively related to florfenicol and thiamphenicol residue in it, after comparing with the Standard Curve, multiplied by the dilution factor, florfenicol and thiamphenicol residue quantity in the sample can be calculated.

# **Reagents And Materials Provided**

- 1. Microtiter plate with 96 wells coated with antigen
- 2. Standard solutions: (6×1ml/bottle),0ppb, 0.05ppb,0.15ppb, 0.45ppb,1.35ppb,4.05ppb
- 3. Spiking standard solution: (1ml/bottle) 1ppm
- 4. Concentrated Enzyme conjugate: 1ml, red cap
- 5. Enzyme conjugate diluent: 10ml, green cap
- 6. Substrate solution A: 7ml, white cap
- 7. Substrate solution B: 7ml, red cap
- 8. Stop solution: 7ml, yellow cap
- 20xconcentrated wash solution: 40ml, transparent cap
- 10. 2xconcentrated extraction solution: 50ml, blue cap
- 11. Extraction buffer: 7ml, red cap

# **Materials Required But Not Supplied**

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



#### Cat: DEIA039

#### 1. Equipments

- Microtiter plate spectrophotometer (450nm/630nm) 1.
- 2. Rotary evaporator or nitrogen drying instruments
- 3. Homogenizer
- 4. Shaker
- 5. Vortex mixer
- 6. Centrifuge
- 7. Analytical balance (inductance: 0.01g)
- 8. Graduated pipette: 10ml
- 9. Rubber pipette bulb
- 10. Volumetric flask: 500ml
- 11. Polystyrene Centrifuge tubes:2ml, 50ml
- 12. Micropipettes: 20µl-200µl, 100µl-1000µl, 250µl-multipipette

#### 2. Reagents

Acetonitrile (AR)

N-hexane (AR)

Trichloroacetic acid (TDA, AR)

Deionized water

# Storage

Storage condition: 2-8°C.

Storage period: 12 months.

# **Specimen Collection And Preparation**

## Notice and precautions before operation:

Please use one-off tips in the process of experiment, and change the tips when absorb different reagent.

Make sure that all experimental tools are clean.

Keep untreated samples in freeze.

Treated samples can't be stored and should be used immediately.

# 2. High detection limit tissue (chicken, duck, pork, chicken liver, pig liver, fish and shrimp), honey, egg

Homogenize the sample with homogenizer;

Weigh 1.0±0.05g homogenized sample into a 50ml polystyrene centrifuge tube, then add 4ml of acetonitrile, vortex for 5min, then centrifuge for separation: 5min / 4000r/min / ambient temperature;

Transfer 1ml of the supernatant organic phase into a 10ml clean glass test tube, dry with 50-60°C water bath under nitrogen gas stream;

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221

Add 1ml of n-hexane, vortex for 30s, then add 1ml of extraction solution(solution 2), vortex for 1min, then transfer the solution into a 2ml polystyrene centrifuge tube, then centrifuge for separation: 5min / 4000r/min / ambient temperature.

Remove the upper layer of organic phase, take 50µl of the lower layer phase for assay;

Dilution factor: 4

## 3. Low detection limit tissue (chicken, duck, pork, chicken liver, pig liver, fish and shrimp), egg

Homogenize the sample with homogenizer;

Weigh 1.0±0.05g homogenized sample into a 50ml polystyrene centrifuge tube, then add 4ml of deionized water, add 50µl of extraction buffer, vortex for 3min, then centrifuge for separation: 5min / 4000r/min / ambient temperature;

Transfer 50µl of supernatant into a 2ml polystyrene centrifuge tube, add 950µl of extraction solution (solution 2), vortex for 30s.

Take 50µl for assay;

Dilution factor: 100

#### 3. Milk (raw milk, UHT pure milk)

Take 1ml of milk sample into a 2ml polystyrene centrifuge tube, then add 1ml of 1% trichloroacetic acid solution(Solution 1), vortex for 2min, then centrifuge for separation: 5min / 4000r/min / ambient temperature;

Transfer 200µl of the supernatant (avoid the floating substance) into a 2ml polystyrene centrifuge tube, add 800µl of extraction solution(Solution 2), mix completely.

Take 50µl for assay;

Dilute factor: 10

#### 4. Feed sample

Weigh 1.0±0.05g of the feed sample into a 50ml polystyrene centrifuge tube, add 10ml of wash solution(Solution 3), vortex for 2min, then centrifuge for separation: 5min / 4000r/min / ambient temperature;

Transfer 50µl of the supernatant(avoid the floating substance) into a 2ml polystyrene centrifuge tube, add 950µl of extraction solution(Solution 2), mix completely.

Take 50µl per well for assay.

**Dilution factor: 200** 

# **Reagent Preparation**

Solution 1: 1% trichloroacetic acid solution

Dissolve 1.0g of trichloroacetic acid with deionized water and to 100ml.

Solution 2: Extraction solution

Dilute the 2x concentrate extraction solution with deionized water in the volume ratio of 1:1 (or according to the requirement), which will be used for sample dilution. The diluted extraction solution can be conserved for 1 month at 4°C.

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

**Fax:** 1-631-938-8221

#### Solution 3: Wash solution

Dilute the 20x concentrated wash solution with deionized water in the volume ratio of 1:19 (or according to the requirement), which will be used to rinse the plate. The diluted wash solution can be conserved for 1 month at 4°C..

## **Assay Procedure**

### Notice before assay

Make sure all reagents and microwells are all at room temperature (20-25°C).

Return all the rest reagents to 2-8°C immediately after used.

Washing the microwells correctly is an important step in the process of assay; it is the vital factor to the reproducibility of the ELISA analysis.

Avoid the light and cover the microwells during incubation.

#### **Assay Steps**

- Take all reagents out at room temperature (20-25°C) for more than 30min, shake gently before use.
- 2. Get the microwells needed out and return the rest into the zip-lock bag at 2-8°C immediately.
- 3. All reagents should be rewarmed before use.
- 4. Number: Number every microwell position and all standards and samples should be run in duplicate. Record the standards and samples positions.
- Dilute the concentrated enzyme conjugate: Dilute the concentrated enzyme conjugate(Kit provided) with the 5. enzyme conjugate diluent(Kit provided) in the volume ratio 1:10(e.g. 0.5ml of concentrated enzyme conjugate + 5ml of enzyme conjugate diluent), mix thoroughly, ready for use. Note: AAA The diluted enzyme conjugate can not be preserved, use immediately.
- Add standard solution/sample and diluted enzyme conjugate: Add 50µl of standard solution(Kit provided) or prepared sample to corresponding wells, add 50µl of diluted enzyme conjugate solution(see step 5) to each well, and then mix gently by shaking the plate manually and incubate for 30min at 25°C with cover.
- Wash: Remove the cover gently and pour the liquid out of the wells and rinse the microwells with 250µl diluted wash solution (solution 3) at interval of 10s for 4-5 times. Absorb the residual water with absorbent paper.
- Coloration: Add 50µl of solution A(Kit provided) and 50µl of solution B(Kit provided) to each well. Mix gently by shaking the plate manually and incubate for 15 min at 25°C with cover(see Precautions).
- Measure: Add 50µl the stop solution(Kit provided) to each well. Mix gently by shaking the plate manually and measure the absorbance at 450nm against an air blank (It's suggested measure with the dualwavelength of 450/630nm. Read the result within 5min after addition of stop solution. We can also measure by sight without stop solution in short of the ELISA reader)

#### Calculation

## 1. Percentage absorbance

The mean values of the absorbance values obtained for the standards and the samples are divided by the absorbance value of the first standard (zero standard) and multiplied by 100%. The zero standard is thus

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

**Fax:** 1-631-938-8221

Cat: DEIA039

made equal to 100% and the absorbance values are quoted in percentages.

Absorbance (%) = 
$$\frac{B}{B_0} \times 100\%$$

B ——absorbance standard (or sample)

B0 ——absorbance zero standard

#### 2. Standard Curve

To draw a standard curve: take the absorbance value of standards as y-axis, semi logarithmic of the concentration of the florfenicol and thiamphenicol standards solution (ppb) as x-axis.

The florfenicol and thiamphenicol concentration of each sample(ppb), which can be read from the calibration curve, is multiplied by the corresponding dilution factor of each sample followed, and the actual concentration of sample is obtained.

#### Please notice:

For data reduction of the ELISA kits, special software has been developed, which can be provided onrequest.

## **Performance Characteristics**

## **Accuracy**

Tissue, honey, egg: 100±30%

Milk: 100±30% Feed: 100±30%

# **Precision**

Variation coefficient of the ELISA kit is less than 10%.

### **Detection Limit**

High detection limit tissue, honey, egg: 0.2ppb

Low detection limit tissue, egg: 5ppb

Milk: 0.5ppb Feed: 10ppb

# Sensitivity

0.05ppb

# **Specificity**

Florfenicol: 100%

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221

Cat: DEIA039

Thiamphenicol: 60% Florfenicol amine: <1% Chloramphenicol: <1%

#### **Precautions**

- The mean values of the absorbance values obtained for the standards and the samples will be reduced if the reagents and samples have not been regulated to room temperature (20-25°C).
- 2. Do not allow microwells to dry between steps to avoid unsuccessful repetitiveness and operate the next step immediately after tap the microwells holder.
- 3. Shake each reagent gently before use.
- 4. Keep your skin away from the stop solution for it is 2M H<sub>2</sub>SO<sub>4</sub>.
- 5. Don't use the kits out of date. Don't exchange the reagents of different batches, or else it will drop the sensitivity.
- Keep the ELISA kits at 2-8°C, do not freeze. Seal rest microwell plates Avoid straight sunlight for the 6. standard sample and the colorless chromogenic reagent are sensitive to light.
- 7. Substrate solution should be abandoned if it turns colors. The reagents may turn bad if the absorbance value (450/630nm) of the zero standard is less than 0.5(A450nm<0.5).
- The coloration reaction needs 10-15min after adding Solution A and Solution B. And you can prolong the incubation time ranges to 20 min or more if the color is too light to be determined, never exceed 30min, On the contrary, shorten the incubation time properly.
- The optimal reaction temperature is 25°C. Higher or lower temperature will lead to the changes of sensitivity and absorbance values

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221