



User's Manual

Human SDMA ELISA Kit

REF

DEIA10342



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human SDMA (Symmetric Dimethylarginine) ELISA Kit is intended for the quantitative determination of SDMA (Symmetric Dimethylarginine) in human serum or plasma. The Human SDMA (Symmetric Dimethylarginine) ELISA Kit is for research use only and not to be used in clinical, therapeutic or diagnostic procedures.

General Description

Dosing of most drugs must be adapted in renal insufficiency, making accurate assessment of renal function an essential component of diagnostics in clinical medicine. Furthermore, even modest impairment of renal function has been recognized as a cardiovascular risk factor. As the most commonly used marker of renal excretory function, serum creatinine concentration, does not adequately respond to mild to moderate impairment of renal function, more sensitive markers for renal excretory function are urgently sought, especially in mild stages of renal impairment. SDMA is a methylated derivative of the amino acid L-arginine (symmetric dimethylarginine). SDMA is eliminated from the body exclusively by renal excretion; therefore SDMA plasma concentration is tightly related to renal function. Thus, quantification of plasma SDMA is an adequate means to assess renal function, as could be demonstrated in a series of recent clinical trials: In 18 clinical studies involving more than 2,100 patients systemic SDMA concentrations were highly correlated with inulin clearance as well as with various clearance estimates and better corresponded to mild renal function impairment than serum creatinine.

Thus, SDMA exhibits properties of a reliable marker of renal function. Furthermore, there is evidence showing that elevated SDMA levels, as they may occur in renal function impairment, may prospectively indicate future risk of cardiovascular disease and mortality independently of the level of renal impairment.

Principles of Testing

The competitive SDMA human ELISA uses the microtiter plate format. SDMA is bound to the solid phase of the microtiter plate. SDMA in the samples is acylated and competes with solid phase bound SDMA for a fixed number of rabbit anti-SDMA antiserum binding sites. When the system is in equilibrium, free antigen and free antigen-antiserum complexes are removed by washing. The antibody bound to the solid phase SDMA is detected by anti-rabbit / peroxidase. The substrate TMB / peroxidase reaction is monitored at 450 nm. The amount of antibody bound to the solid phase SDMA is inversely proportional to the SDMA concentration of the sample.

Reagents And Materials Provided

1. MT-Strips (12 strips): 8 wells each, break apart, precoated with SDMA
2. Standards 1-6 (6 vials): each 4 ml, ready for use, concentration as below:

Standard	1	2	3	4	5	6
μmol/l	0	0.2	0.4	0.7	1.2	3.0
Ng/ml	0	40	81	141	242	606

3. Control 1&2 (2 vials): each 4 ml, ready for use. Range: see q.c; certificate
4. Acylation Buffer (1 vial): 3.5 ml, ready for use, blue coloured
5. Acylation Reagent (3 vials): lyophilized, dissolve contents. In 3 ml solvent before use
6. Antiserum (1 vial): 7 ml, ready for use, yellow coloured. Rabbit-anti-N-acyl-SDMA
7. Enzyme Conjugate (1 vial): 13 ml, ready for use. Goat anti-rabbit-IgG-peroxidase
8. Wash Buffer (1 vial): 20 ml, concentrated (50x). Dilute contents with dist. water to 1000 ml total volume
9. Substrate (1 vial): 13 ml TMB solution, ready for use
10. Stop Solution (1 vial): 13 ml, ready for use. Contains 0.3 M sulphuric acid, not corrosive
11. Reaction Plate (1 piece): for acylation
12. Equalizing Reagent (1 vial): lyophilized, dissolve contents with 21 ml dist. water. Dissolve carefully to minimize foam formation
13. Solvent (2 vials): 5 ml, contains DMSO. Please note that solvent reacts with many plastic materials including plastic trays. Solvent does not react with normal pipette tips and with glass devices
14. Foil (2 pieces)

Materials Required But Not Supplied

1. Pipettes 20, 50, 100 and 200 µl
2. Multipipette
3. Orbital shaker
4. Microplate washing device
5. Microplate photometer (450 nm)
6. Vortex mixer, roll mixer

Storage

On arrival, store the kit at 2-8 °C. Once opened the kit is stable until its expiry date.

For stability of prepared reagents refer to Preparation of Reagents.

Do not use components beyond the expiration date shown on the labels.

Do not mix various lots of any kit component within an individual assay.

Specimen Collection And Preparation

Serum and Plasma

The test can be performed with serum as well as with EDTA plasma.

Hemolytic and lipemic samples should not be used.

The samples can be stored up to 6 hours at 2 - 8 °C. For a longer storage (up to 18 months) the samples must be kept frozen at -20 °C

Repeated freezing and thawing should be avoided.

PREPARATION OF SAMPLES (Acylation)

The wells of the reaction plate for the acylation can be used only once. Please mark the respective wells before use to avoid repeated use.

1. Pipette each 20 µl standard 1 - 6, each 20 µl control 1 & 2 and each 20 µl patient sample into the respective wells of the Reaction Plate.
2. Pipette 20 µl Acylation Buffer into all wells.
3. Pipette 200 µl reconstituted Equalizing Reagent into all wells.
4. Mix the reaction plate for 10 seconds.
5. Prepare Acylation Reagent freshly and pipette 50 µl prepared Acylation Reagent each into all wells, mix immediately. It is recommended to use a multipipette, fill it directly from the vial and add the Acylation Reagent to the wells. Colour changes to violet.
6. Incubate for 20 minutes at room temperature (approx. 20 °C) on an orbital shaker. Do not cover wells or plate, leave the plate open on the shaker.

Take each 25 µl of the acylated samples for the SDMA human ELISA.

Reagent Preparation

Microtiter strips

Before opening the packet of strip wells, allow it to stand at room temperature for at least 10 minutes. After opening, keep any unused wells in the original foil packet with the desiccant provided. Reseal carefully and store at 2-8 °C.

Wash Buffer

Dilute the contents with dist. water to a total volume of 1000 ml, mix shortly. The diluted wash buffer must be stored at 2 - 8 °C and is stable for 4 weeks. For longer storage the diluted wash buffer has to be stored frozen at -20°C.

Equalizing Reagent

Dissolve the contents with 21 ml dist. water, mix shortly and leave on a roll mixer for 20 minutes. Avoid excess formation of foam. The reconstituted Equalizing Reagent should be stored frozen at -20 °C and is stable until expiry date.

Acylation Reagent

Dissolve the contents of one bottle in 3 ml Solvent and shake for 10 minutes on an orbital shaker. After use the reagent has to be discarded. The Acylation Reagent has always to be prepared immediately before use and is stable for minimum 3 hours. The two other bottles allow a second and third run of the test. If the whole kit is to be used in one run it is recommended to pool the dissolved contents of the two vials of Acylation Reagent.

Please note that Solvent reacts with many plastic materials including plastic trays which are used as reservoir for multichannel pipettes. Solvent does not react with normal pipette tips and with glass devices. It is recommended to use a multipipette, fill it directly from the vial and add the Acylation Reagent to the wells.

Assay Procedure

Bring all reagents to room temperature and mix them carefully, avoid development of foam.

1. Sample Incubation

Pipette each 25 µl prepared Standards 1 to 6, 20 µl prepared controls and 25 µl prepared samples into the respective wells of the coated microtiter strips (duplicates are recommended).

Pipette each 50 µl Antiserum into all wells.

Cover the plate with adhesive foil and incubate Microtiter Strips for 90 minutes at room temperature (20 - 25°C) on an orbital shaker.

2. Washing

Discard or aspirate the contents of the wells and wash thoroughly with each 300 µl Wash Buffer. Repeat the washing procedure 4 times. Remove residual liquid by tapping the inverted plate on clean absorbent paper.

3. Conjugate Incubation

Pipette each 100 µl enzyme conjugate into all wells.

Incubate for 30 minutes at room temperature on an orbital shaker.

4. Washing Repeat step 2.

5. Substrate Incubation

Pipette each 100 µl Substrate into all wells and incubate for 25 ± 5 minutes at room temperature on an orbital shaker.

6. Stopping

Pipette each 100 µl Stop Solution into all wells.

7. Reading


Read the optical density at 450 nm (reference wavelength between 570 and 650 nm) in a microplate photometer.

Calculation

On a semilogarithmic graph paper the concentration of the standards (x-axis, logarithmic) are plotted against their corresponding optical density (y-axis, linear). Cubic spline, 4 parameter or similar iteration procedures are recommended for evaluation of the standard curve.

The concentration of the controls and samples can be read directly from this standard curve by using their average optical density.

Typical Standard Curve



Conversion factor: 1 µmol/l = 202 ng/ml = 20.2 µg/dl

Reference Values

Serum, EDTA-Plasma

Humans : 0.30 – 0.75 µmol/l (6.0 – 15 µg/dl)

The reference ranges given above should only be taken as a guideline. It is recommended that each laboratory should establish its own reference values.

Precision

	Range (µmol/l)	Intra Assay CV
EDTA-Plasma	0.52-0.82	6.2-4.9 %

	Range (µmol/l)	Inter Assay CV
EDTA-Plasma	0.52-1.21	2.0-8.8 %

Sensitivity

0.03 µmol/l

Specificity

Substance	Cross Reactivity (%)
SDMA	100
ADMA	0.74
NMMA	0.76
Homoarginine	0.04
Arginine	0.01

Linearity

	Range (µmol/l)	Highest Dil.	Mean %	Range %
EDTA-Plasma	0.23-1.72	1:6 with water	97	89-105

Recovery

	Range (µmol/l)	Mean %	Range %
EDTA-Plasma	0.43-1.44	97	86-104
Serum	0.45-1.72	93	88-102

Precautions

1. For research use only
2. Disposable gloves should be used.
3. Material of animal origin used in the preparation of the kit has been obtained from animals certified as healthy, but these materials should be handled as potentially infectious.