



User's Manual

Chikungunya Virus IgG ELISA Kit



DEIA2162



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The ELISA Chikungunya Virus IgG test is qualitative and quantitative immunoassay for the detection of human antibodies in serum or plasma directed against Chikungunya Virus.

General Description

Chikungunya Virus (CHIKV) is a member of the Alphavirus genus. It is also classified as an Arbovirus due to vector-dependent transmission by mosquitos. The enveloped virus carries structural glycoproteins E1-E3 which are immunogenic. Several non-structural proteins (nsP) and a ss-(+) RNA genome are present within the viral capsid.

Chikungunya fever is characterized by fever, joint pain and swelling, rash, and flu-like symptoms. Symptoms usually appear within a period of 1-2 weeks post-infection. The disease is characterized by a biphasic course with an initial stage of high virus load and a convalescent stage without viremia. In rare cases, neurological disorders like Guillain-Barré syndrome and meningoencephalitis have occurred.

Principles of Testing

The ELISA (Enzyme-Linked Immunosorbent Assay) is an immunoassay suited to the detection of antibodies. The reaction is based on the specific interaction of antibodies with their corresponding antigen. The test strips of the ELISA microtiter plate are coated with specific antigens of the pathogen of interest. If antibody in a sample is present, they bind to the fixed antigen. A secondary antibody, which has been conjugated with the enzyme alkaline phosphatase, detects and binds to the antigen-antibody complex. The colorless substrate pnitrophenylphosphate is then converted into the colored product p-nitrophenol. The signal intensity of this reaction product is proportional to the concentration of the antibody in the sample and is measured photometrically.

Reagents And Materials Provided

- 1. Break apart microtiter test strips each with eight antigen coated single wells**, (altogether 96), 1 frame. The coating material is inactivated. 12 pieces
- 2. Standard serum (ready-to-use)**, Human serum in protein-containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface Antigen) and anti-HCV Ab; Preservative: <0.1% sodium azide; coloring: Amaranth O. 2 x 2 ml
- 3. Negative control serum (ready-to-use)**, Human serum in protein-containing phosphate buffer; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface Antigen) and anti-HCV Ab; Preservative: <0.1% sodium azide; coloring: Lissamin Green V. 2 ml
- 4. Anti-human IgG conjugate (ready-to-use)**, Anti-human IgG polyclonal antibody, Conjugated to alkaline phosphatase, stabilized with protein stabilization solution; Preservative: <0.1% methylisothiazolone, <0.1% bromnitrodioxane. 13ml
- 5. Washing solution concentrate (sufficient for 1000ml)**, Sodium chloride solution with Tween 20 and

30mM Tris-HCl, pH 7.4; Preservative: <0.1% sodium azide. 33.3ml

6. Dilution buffer (ready-to-use), Protein-containing phosphate buffer with Tween 20; Preservative: <0.1% sodium azide; coloring: 0.01g/l Bromphenol blue. 2 × 50ml

7. Stopping solution (ready-to-use), <0.1N sodium hydroxide, 40mM EDTA . 15ml

8. Substrate (ready-to-use), Para-nitrophenylphosphate in solvent-free buffer; Preservative: <0.1% sodium azide. 13ml

Materials Required But Not Supplied

1. Common laboratory equipment
2. Photometer for microtiter plates with filter, wavelength 405nm, recommended reference wavelength 620nm-690nm (e.g., 650nm)
3. Microtiter plate washer
4. Incubator 37°C
5. Moist chamber
6. Distilled water
7. Optional: ELISA control

Storage

Reagent	Storage	Stability
Microtiter strips (coated with antigen)	Unopened	See expiry date
	After opening at 2-8°C in closed aluminum bag with desiccant	6 months
Control sera / Standard sera	Unopened	See expiry date
	After opening at 2-8°C	6 months
Conjugate	Unopened	See expiry date
	After opening at 2-8°C	6 months
Dilution buffer	Unopened	See expiry date
	after opening at 2-8°C	6 months
Washing solution	Unopened / after opening at 2-8°C	See expiry date
	Working dilution at 2-8°C	2 weeks
	Working dilution at room temperature	1 week
Substrate	Unopened	See expiry date
	After opening at 2-8°C	6 months
Stopping solution	Unopened	See expiry date
	After opening at 2-8°C	6 months

Specimen Collection And Preparation

1. Sample Preparation and Storage

Lipaemic, hemolytic or icteric samples (serum or plasma) should only be tested with caution. Obviously contaminated samples should not be tested. Serum or plasma (EDTA, citrate, heparin) collected according to

standard laboratory methods are suitable samples. Samples must not be thermally inactivated.

Dilution of Samples

Before running the test, samples (V_1) must be diluted in dilution buffer (V_2) as follows:

$V_1 + V_2 = 1:500$	add	10ul	sample
	each to	1000ul	dilution buffer (= 1:100)
	each to	50ul	of the first dilution
		200ul	dilution buffer (= 1:5)

After dilution and before pipetting into the microtiter plate, the samples must be mixed thoroughly to prepare a homogenous solution.

2. Sample Storage

Samples should not be stored for more than 7 days at 2-8°C. Extended storage is possible at $\leq -20^\circ\text{C}$. Avoid repeated freezing and thawing of samples. Diluted samples can be stored at 2-8°C for one week.

Reagent Preparation

Bring all reagents to room temperature before testing.

1. Microtiter Test Strips

The microtiter test strips labeled with abbreviations for pathogen and immunoglobulin class are packed with a desiccant in an aluminum bag. To open the aluminum bag of the microtiter plate, please cut off the top of the marked side only in order to guarantee proper resealing. Take unrequired wells out of the frame and put them back into the aluminum bag. Close bag carefully to ensure airtight conditions. Do not use strips if the aluminum bag is damaged or if the bag with remaining strips and desiccant was not properly resealed.

2. Negative Control Sera / Standard Sera (ready-to-use)

Negative control and standard sera are ready-to-use and must not be diluted any further. For each test run (independent of the number of microtiter test strips to be used) negative control and standard sera must be included. Standard sera should be set up in duplicate. Do not treat negative control and standard sera with Rf-absorbent.

3. Anti-human IgG AP-Conjugate (ready-to-use)

The required conjugate concentration (i.e., +, ++, +++) is indicated on the quality control certificate. Please refer also to the specification on the label. Avoid contamination.

4. Washing Solution (Concentrate)

Dilute washing buffer concentrate (V_1) 1:30 with distilled H_2O to a final volume of V_2 . Bottles used for the working dilution should be cleaned regularly. Discard cloudy solutions.

Example:

Buffer concentrate (V_1)	Final volume (V_2)
33.3ml	1000ml
1.0ml	30ml

5. Dilution Buffer for Samples (ready-to-use)

Discard cloudy solutions.

6. Substrate (ready-to-use)

Substrate in unopened bottle may have a slight yellow color which does not reduce the quality of the product! Avoid contamination.

7. Stopping Solution (ready-to-use)

Assay Procedure

1. Place the required number of wells in the frame and prepare a protocol sheet.
2. Add each 100ul of diluted sample or ready-to-use negative control / standard sera into the appropriate wells of microtiter test strips. Spare one well for substrate blank, e.g.:

Well	Quantitative ELISA
A1	Substrate blank
B1	Negative control
C1	Standard serum
D1	Standard serum
E1	Sample 1 . . .
F1	Sample 2 . . .

3. Sample incubation for 60 minutes (+/- 5 min) at 37°C (+/- 1°C) in moist chamber.
4. After incubation wash all wells with washing solution (by automated washer or manually):
aspirate or shake out the incubation solution
fill each well with 300ul washing solution
aspirate or shake out the washing solution
repeat the washing procedure 3 times (altogether 4 times!)
dry by tapping the microtiter plate on a paper towel
5. Addition of conjugate: Add 100ul of the ready-to-use IgG conjugate to the appropriate wells (except substrate blank).
6. Conjugate incubation for 30 minutes (+/- 1 min) at 37°C (+/- 1°C) in moist chamber.
7. After incubation wash all wells with washing solution (see above).
8. Addition of substrate: Add 100ul of ready-to-use substrate solution to each well (including well for substrate blank!)
9. Substrate incubation for 30 minutes (+/- 1 min) at 37°C (+/- 1°C) in moist chamber. Ensure incubation is in the dark.
10. Stopping the reaction: Add 100ul of stopping solution to each well, shake microtiter plate gently to mix.
11. Read extinction: Read optical density (OD) within 60 minutes at 405nm against substrate blank, reference



wavelength between 620nm and 690nm (e.g. 650nm).

Calculation

1. ELISA Chikungunya Virus IgG

The mathematical curve fitting for antibody quantification with ELISA immunoassays is based on the 4-parameter logistic (4 PL) function.

$$Activity (U/ml) = e^{C - \frac{1}{B} \ln\left(\frac{D-A}{OD(Patient) * F - A}\right)}$$

The 4 parameters A, B, C, and D are representative for the exact shape of the standard curve:

Parameter A: Lower asymptote (OD)

Parameter B: Slope of the curve

Parameter C: Inflection point

Parameter D: Upper asymptote (OD)

CD establishes a lot-specific 4 PL standard curve for this ELISA immunoassay in multiple test runs under optimal test conditions. The four parameters are indicated on the quality control certificate of each individual ELISA test.

For the adaptation of the test level to the given 4 PL standard curve, the correction factor F is calculated by dividing the standard reference OD value indicated on the quality control certificate with the measured, and consequently test run-specific, standard OD value.

$$F = \frac{\text{STD reference OD value}}{\text{measured STD OD value}}$$

By multiplying the OD values obtained from samples with the correction factor F, the level of each individual test run is adjusted to the given 4 PL standard curve. Thereby, interassay deviations are compensated for and antibody activities can be directly evaluated from the 4 PL standard curve.

After subtraction of the substrate blank from all measured OD values and calculation of the mean OD value of the standard serum (STD), tested in duplicate, the evaluation of antibody activities from the optical measurement signals (OD) of samples can be performed with 4PL function presented above.

2. Borderline Ranges

The borderline ranges of the ELISA Chikungunya Virus IgG test are specified on the quality control certificates and indicate the range of borderline test results.

Values below this range indicate a negative result; values above the borderline range indicate a positive result.

3. Limits of Quantification

The limits of quantification are specified on the quality control certificate of the ELISA Chikungunya Virus IgG. The linearity of dilution within this range has been demonstrated in comprehensive evaluation studies. If a

sample shows a test result above the upper limit of quantification, the sample may be tested at a higher dilution. The resulting antibody activity must then be multiplied by the additional dilution factor.

4. Qualitative Evaluation with ELISA Chikungunya Virus IgG

For the ELISA test evaluation, a lot-specific quality control certificate with standard curve and an evaluation table is included in the test kit so that the obtained OD values may be assigned to the corresponding antibody activities. The substrate blank must be subtracted from all OD values prior to evaluation. Mean OD value of the standard serum (STD), tested in duplicate, has to be used.

Method 1:

In the first line of the evaluation table, several ranges of OD values for the standard serum that cover the entire standard validity range are provided. Choose the column that corresponds to the measured mean OD value obtained for the standard serum. This column contains the upper and lower cut-off OD values for evaluation of antibody in the sample. OD values below the lower cut-off are considered negative values, and OD values above the upper cut-off are considered positive values. Use of the correction factor F is not needed for this evaluation method.

Method 2:

To fix the cut-off ranges, multiply the mean value of the OD measured for the standard serum with the numerical data of the quality control certificate (see special case formulas), e.g.:

$OD = 0.502 \times MW(STD)$ with upper cut-off

$OD = 0.352 \times MW(STD)$ with lower cut-off

Calculation example: Standard serum mean OD = 0.64

Upper cut-off: $OD = 0.502 \times 0.64 = 0.321$

Lower cut-off: $OD = 0.352 \times 0.64 = 0.225$

If the measured OD of the standard serum is 0.64, the range of the cut-off is between 0.225-0.321.

5. Criteria of Validity

- The substrate blank must be <0.25 OD.
- The negative control must be negative.
- The mean OD value (after subtraction of the substrate blank!) of the standard serum must be within the validity range which is given on the lot-specific quality control certificate.
- The variation of OD values of the standard serum or cut-off serum must not be higher than 20%. If these criteria are not met, the test is not valid and must be repeated.

Specificity

Potential cross-reactivities of the ELISA Chikungunya Virus IgG were evaluated with known positive sera (10 each) for Zika Virus IgG, Dengue Virus IgG, Rubella Virus IgG, TBE Virus IgG, Influenza A Virus IgG, Epstein-Barr Virus VCA IgG, West Nile Virus IgG, Parvovirus B19 IgG, Measles Virus IgG, Mumps Virus IgG, Varicella-Zoster Virus IgG, and Influenza B Virus IgG as well as sera positive for rheumatoid factor (RF) and anti-nuclear antibodies (ANA). Cross-reactivity with ten Zika Virus IgG sera were observed. Other crossreactivities cannot be ruled out in general.

Interferences

No interferences have been detected for sera with concentrations up to 2.00g/L hemoglobin, 11.50g/L lipemia/triglyceride, or 0.201g/L bilirubin (conjugated and unconjugated).

Precautions

1. Evidence of Deterioration

Optimum results can only be achieved if the instructions are strictly followed. Only use ELISA reagents when using ELISA immunoassays. The components must not be exchanged for reagents of other manufacturers. Standard and control sera of ELISA immunoassays are defined exclusively for the test kit to be used and must not be used in other lots.

Each ELISA classic test contains a ready-to-use sample dilution buffer. In some cases the use of special dilution buffers is necessary to guarantee consistent quality and reliable results. There are three different conjugate concentrations for each immunoglobulin class (IgA, IgG, IgM) indicated on the label as + (low), ++ (medium), and +++ (high). Dilution or alteration of the reagents may result in a loss of sensitivity. Use aseptic techniques when removing aliquots from the reagent tubes to avoid contamination.

Reproducibility of test results is dependent on thorough mixing of the reagents. Agitate the vials containing control sera before use and also all samples after dilution (e.g., by using a vortex mixer).

Be sure to pipette carefully and comply with the given incubation times and temperatures.

Significant time differences between pipetting the first and last well of the microtiter plate when dispensing samples and control sera, conjugate or substrate can result in different preincubation times, which may influence the precision and reproducibility of the results. Avoid exposure of reagents to strong light during storage and incubation.

Adequate washing avoids test unspecificities. Therefore, the washing procedure should be carried out carefully. All of the flat bottom wells should be filled with equal volumes of washing buffer. At the end of the procedure ensure that the wells are free of all washing buffer in order to avoid uncontrolled dilution effects. Avoid foaming!

Reagents must be tightly closed after use to avoid evaporation and contamination. Take care not to mix up the caps of the bottles and/or vials.

The ELISA immunoassay is only valid if the lot-specific validation criteria on the quality control certificate are fulfilled.

Statements of Warning

The ELISA is designed for use by qualified personnel who are familiar with good laboratory practice. All kit reagents and human samples should be handled carefully using established good laboratory practice.

1. This kit contains human blood components. Although all control- and cut-off sera have been tested and found negative for anti-HIV Ab, HBs-Ag (Hepatitis B Virus surface Antigen) and anti-HCV Ab, they should be considered potentially infectious.
2. Do not pipette by mouth.
3. Do not smoke, eat, or drink in areas in which samples or kit reagents are handled.
4. Wear disposable gloves, laboratory coat, and safety glasses while handling kit reagents or samples. Wash

hands thoroughly afterwards.

5. Samples and other potentially infectious material should be decontaminated after use.
6. Reagents should be stored safely and be inaccessible to unauthorized access, e.g. children.