



One Step Malaria Pf/Pan (Whole blood) Test (IVDTS003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use

One Step Malaria Pf/Pan Test is a lateral flow chromatographic immunoassay for the detection of Malaria plasmodium lactate dehydrogenase (pLDH) in vitro assay. It is used for detection of Malaria P. falciparum and other Plasmodium (P.vivax, P. malariae, P. ovale) in human blood. This device is intended to be used as a screening test and as an aid in the diagnosis of infection with Plasmodium. Any reactive specimen with One Step Malaria Pf/Pan Test must be confirmed with alternative testing method(s) and clinical findings.

Storage

The test kit can be stored at 4 to 30°C in the sealed pouch to the date of expiration. The test kit should be kept away from direct sunlight, moisture and heat.

45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 Fax: 1-631-938-8221

Email:info@creative-diagnostics.com

© Creative Diagnostics All Rights Reserved

1/1