



**User's Manual** 

# Human GRN ELISA Matched Antibody Pair



**ABPR-0394** 





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

#### **Creative Diagnostics**

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221 Email: info@creative-diagnostics.com Web: www.creative-diagnostics.com

## PRODUCT INFORMATION

#### **Intended Use**

This antibody pair set comes with matched antibody pair to detect and quantify protein level of human GRN.

## **General Description**

Granulins are a family of secreted, glycosylated peptides that are cleaved from a single precursor protein with 7.5 repeats of a highly conserved 12-cysteine granulin/epithelin motif. The 88 kDa precursor protein, progranulin, is also called proepithelin and PC cell-derived growth factor. Cleavage of the signal peptide produces mature granulin which can be further cleaved into a variety of active, 6 kDa peptides. These smaller cleavage products are named granulin A, granulin B, granulin C, etc. Epithelins 1 and 2 are synonymous with granulins A and B, respectively. Both the peptides and intact granulin protein regulate cell growth. However, different members of the granulin protein family may act as inhibitors, stimulators, or have dual actions on cell growth. Granulin family members are important in normal development, wound healing, and tumorigenesis.

### Reagents And Materials Provided

Antibody pair set content:

- Capture antibody: mouse monoclonal anti-GRN (100 µg)
- Detection antibody: rabbit purified polyclonal anti-GRN (50 μg) 2.

Note: Reagents are sufficient for at least 3-5 x 96 well plates using recommended protocols.

## **Reconstitution And Storage**

Store reagents of the antibody pair set at -20°C or lower. Please aliquot to avoid repeated freeze thaw cycle. Reagents should be returned to -20°C storage immediately after use.

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

**Fax:** 1-631-938-8221

Email: info@creative-diagnostics.com