



# Human HLA-DPA1 ELISA Matched Antibody Pair



**ABPR-0425** 





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

### **Creative Diagnostics**

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221 Email: info@creative-diagnostics.com Web: www.creative-diagnostics.com

## PRODUCT INFORMATION

#### **Intended Use**

This antibody pair set comes with matched antibody pair to detect and quantify protein level of human HLA-DPA1.

## **General Description**

HLA-DPA1 belongs to the HLA class II alpha chain paralogues. This class II molecule is a heterodimer consisting of an alpha (DPA) and a beta (DPB) chain, both anchored in the membrane. It plays a central role in the immune system by presenting peptides derived from extracellular proteins. Class II molecules are expressed in antigen presenting cells (APC: B lymphocytes, dendritic cells, macrophages). The alpha chain is approximately 33-35 kDa and its gene contains 5 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the two extracellular domains, exon 4 encodes the transmembrane domain and the cytoplasmic tail. Within the DP molecule both the alpha chain and the beta chain contain the polymorphisms specifying the peptide binding specificities, resulting in up to 4 different molecules.

# Reagents And Materials Provided

Antibody pair set content:

- Capture antibody: rabbit MaxPab® affinity purified polyclonal anti-HLA-DPA1 (100 ug)
- Detection antibody: mouse purified polyclonal anti-HLA-DPA1 (20 ug)

Note: Reagents are sufficient for at least 1-2 x 96 well plates using recommended protocols.

## Reconstitution And Storage

Store reagents of the antibody pair set at -20°C or lower. Please aliquot to avoid repeated freeze thaw cycle. Reagents should be returned to -20°C storage immediately after use.

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

**Fax:** 1-631-938-8221

Email: info@creative-diagnostics.com