



# **Human PSMD4 ELISA Matched Antibody Pair**



**ABPR-0735** 





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

### **Creative Diagnostics**

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221 Email: info@creative-diagnostics.com Web: www.creative-diagnostics.com

## PRODUCT INFORMATION

#### **Intended Use**

This antibody pair set comes with matched antibody pair to detect and quantify protein level of human PSMD4.

# **General Description**

The 26S proteasome is a multicatalytic proteinase complex with a highly ordered structure composed of 2 complexes, a 20S core and a 19S regulator. The 20S core is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. The 19S regulator is composed of a base, which contains 6 ATPase subunits and 2 non-ATPase subunits, and a lid, which contains up to 10 non-ATPase subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. An essential function of a modified proteasome, the immunoproteasome, is the processing of class I MHC peptides. This gene encodes one of the non-ATPase subunits of the 19S regulator lid. Pseudogenes have been identified on chromosomes 10 and 21.

## Reagents And Materials Provided

Antibody pair set content:

- Capture antibody: Magic™ rabbitamp; reg; affinity purified polyclonal anti-PSMD4 (100 µg)
- 2. Detection antibody: mouse monoclonal anti-PSMD4, IgG2a Kappa (20 μg)

Note: Reagents are sufficient for at least 3-5 x 96 well plates using recommended protocols.

## **Reconstitution And Storage**

Store reagents of the antibody pair set at -20°C or lower. Please aliquot to avoid repeated freeze thaw cycle. Reagents should be returned to -20°C storage immediately after use.

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

**Fax:** 1-631-938-8221

Email: info@creative-diagnostics.com