



User's Manual

Human interferon-alpha/IFNA5 Antibody Pair Set



ABPR-ZB247



5 Plates, 15 Plates



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221

Email: info@creative-diagnostics.com

Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Quantitative determination of Human interferon-alpha/IFNA5

General Description

ABPR-ZB247 is a solid phase sandwich ELISA for quantitative determination of Human interferonalpha/IFNA5.

Reagents And Materials Provided

- Capture Antibody: 1 mg/mL of rabbit anti-Human IFNA5/IFNaG monoclonal antibody [CABT-ZB668] (in PBS, pH 7.4).
- Detection Antibody: 0.2 mg/mL of rabbit anti-Human IFNA5/IFNaG polyclonal antibody conjugated to HRP (in PBS, 50 % HRP-Protector, pH 7.4, store at 4°C).
- Standard: Each vial contains 90 ng of recombinant Human IFNA5/IFNaG. 3.

Reconstitution And Storage

Reconstitution

- Capture Antibody: Dilute to a working concentration of 2 µg/mL in PBS before coating.
- 2. Detection Antibbody: Dilute to working concentration of 0.5 μg/mL in detection antibody dilution buffer before
- 3. Standard: Reconstitute with 1 mL detection antibody dilution buffer.

Storage

- Capture Antibody: Aliquot and store at -20°C to -80°C for up to 6 months from date of receipt. Avoid repeated freeze-thaw cycles.
- Detection Antibody: Store at 4°C and protect it from prolonged exposure to light for up to 6 months from date of receipt. DO NOT FREEZE!
- Standard: Store lyophilized standard at -20°C to -80°C for up to 6 months from date of receipt. Aliquot and store the reconstituted Standard at -80°C for up to 1 month. Avoid repeated freeze-thaw cycles.

Detection Range

Assay range: 39.06-2500 pg/mL

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

Fax: 1-631-938-8221

Email: info@creative-diagnostics.com