



Benzodiazepine [HRP] (DAG-WT186)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Benzodiazepine conjugated with HRP
Species	N/A
Conjugate	HRP
Applications	ELISA
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.5 mL
Buffer	0.05M Tris with protein stabilizers
Preservative	0.005% Thimerosal
Storage	Store at 2-8°C, do not freeze.

BACKGROUND

Introduction

Benzodiazepine is the most widely prescribed anxiolytic drug and is used extensively as an anti-anxiety agent, a hypnotic, a muscle relaxant and an anti-convulsant. Benzodiazepine is usually taken orally, or by intravenous or intramuscular injection. High dose abuse of benzodiazepine can result in drowsiness and confusion, and potentiates alcohol and other CNS depressant drugs. Benzodiazepine is metabolized in the liver and excreted in urine. Benzodiazepine and its metabolites are the primary urinary markers for detecting benzodiazepine use. The Benzodiazepines Screen Tests are a rapid, visual, lateral flow, competitive immunochromatographic assay for the qualitative detection of benzodiazepine and its metabolites in human urine. These in vitro diagnostic screening tests are based

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221 © Creative Diagnostics All Rights Reserved

immunoassay principles designed specifically for the assay and identification of benzodiazepine and its metabolites at the cut-off level of 300 ng/ml or higher, as set by SAMHSA.

Keywords

Benzodiazepines; BZO; benzodiazepine; benzo; BZD

Tel: 1-631-624-4882 Fax: 1-631-938-8221

2/2