



Native Bone Alkaline Phosphatase (DAG-WT1923)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	WB, ELISA
Format	Lyophilized
Size	5 µg
Storage	Store at -20°C

BACKGROUND

Introduction	Bone Specific Alkaline Phosphatase is used to diagnose and assess the severity of metabolic bone disease including Paget disease, osteomalacia, and other states of high bone turnover. The test can monitor efficacy of antiresorptive therapies including postmenopausal osteoporosis treatment.
Keywords	Bone Specific Alkaline Phosphatase; BALP