



NO-Tyramine [G-BSA] (DAG3386)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	NO-Tyramine, G-BSA-conjugated
Specificity	NO-Tyramine conjugated with glutaraldehyde (G) and bovine serum albumin (BSA), NO = nitrosylated.
Species	N/A
Conjugate	G-BSA
Applications	immunohistochemistry and immunocytochemistry
Reconstitution	Reconstituted in deionized water (250 μg)
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Keywords Tyramine

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221