



Recombinant Norovirus GII.17 P Domain (DAGA-3045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	Multi-step procedure including affinity chromatography.
Conjugate	Unconjugated
Applications	Suitable for use as a control in Lateral Flow Assays and as an immunogen. Suitable for spiking into stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 mg
Buffer	Tris Buffer
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for < 1 week. For long term storage, aliquot and store at -20°C to avoid multiple freeze/thaw cycles.