



Transferrin protein [FITC] (DAGA-3175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Transferrin protein (FITC)
Purity	Transferrin protein (FITC) was purified by delipidation, selective precipitation and tandem molecular sieve chromatography followed by dialysis.
Conjugate	FITC
Applications	IF
Format	Lyophilized from 0.02M K3PO4, pH 7.2, with 0.15M NaCl, 10 mg/ml BSA and 0.01% NaN3. Immunoglobulin and protease free.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221