



Turkey IgG [FITC] (DAGA-3180)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Turkey Immunoglobulin G (FITC)
Purity	Turkey IgG (FITC) was purified by delipidation, salt fractionation and ion exchange chromatography followed by dialysis.
Conjugate	FITC
Applications	IF
Format	Supplied liquid in 0.01 M Sodium Phosphate, 0.14 M Sodium Chloride, pH 7.6 with 10 mg/ml BSA added as a stabilizer.
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Storage	Store at 4 centigrade until reconstitution. Following reconstitution aliquot and freeze at -20 centigrade for long term storage. Avoid repeated freeze/thaw cycles.