



User's Manual

Porcine Total Nitric Oxide ELISA Kit

REF

DEIA-BJ2834



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Porcine Total Nitric Oxide ELISA Kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Total Nitric Oxide. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

Detection Range

5-100 µmol/L

Sensitivity

1.0 µmol/L