



**User's Manual**

# Human Angiotensin II Receptor 2 ELISA kit



**DEIA-BJ982**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

Human Angiotensin II Receptor 2 ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of the Angiotensin II Receptor 2. This ELISA kit is for research use only, not for therapeutic or diagnostic applications.

### Storage

All components of this kit are stable at 2-8°C until the kit's expiration date.

### Detection Range

0.5-10 ng/mL

### Sensitivity

0.1 ng/mL