



## User's Manual

# AA(Arachidonic Acid) ELISA Kit



DEIA-NS2307-40



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION


### Intended Use

AA ELISA Kit allows for the in vitro quantitative determination of AA concentrations in serum, plasma, tissue homogenates and other biological fluids.

### Storage

2-8 °C for 6 months

### Typical Standard Curve



### Precision

Intra-Assay: CV<8%

Inter-Assay: CV<10%

### Detection Range

1.563-100ug/ml

### Sensitivity

< 0.938ug/ml

### Linearity

The linearity of the kit was assayed by testing samples spiked with appropriate concentration of AA and their serial dilutions. The results were demonstrated by the percentage of calculated concentration to the

expected.

Sample	1:2	1:4	1:8
serum(n=5)	86-103%	94-105%	90-100%
EDTA plasma(n=5)	89-101%	84-96%	91-101%
heparin plasma(n=5)	83-90%	80-100%	82-96%

## Recovery

Matrices listed below were spiked with certain level of AA and the recovery rates were calculated by comparing the measured value to the expected amount of AA in samples.

Matrix	Recovery range(%)	Average(%)
serum(n=5)	85-104	93
EDTA plasma(n=5)	95-103	99
heparin plasma(n=5)	86-97	91