



User's Manual

MGO ELISA kit





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221

PRODUCT INFORMATION

Intended Use

The competitive MGO ELISA kit is to be used for in vitro quantitative determination of MGO in plasma and faeces samples.

Storage

Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.

Detection Range

15.6 - 1000 AU/mL

Detection Limit

15.6 AU/mL

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221

Email: info@creative-diagnostics.com