



User's Manual

Florfenicol & Thiamphenicol ELISA Kit



DEIA-XY2267



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Florfenicol & Thiamphenicol ELISA Kit is intended for the quantitative determination of Thiamphenicol in raw milk, egg, seafood, tissue.

Reagents And Materials Provided

1. Microtiter well plate: 8 microwells/strip, 12 strips/plate.
2. Florfenicol Standards 0, 0.05, 0.1, 0.3, 1, 3 ppb, 1.5 mL each. (Ready to use)
3. 10X Concentrate Extraction/Washing Solution, 50 mL/bottle/kit.
4. Enzyme Conjugate Diluent, 50 mL/bottle/kit.
5. 100X Concentrate Enzyme Conjugate, 100 μ L/bottle/kit.
6. Substrate Solution (TMB), 12 mL/bottle/kit.
7. Stop Solution, 13 mL/bottle/kit.

Detection Limit

The lowest detection limit of the CD Florfenicol & Thiamphenicol ELISA kit is 0.05 ppb. The OD value of this concentration is significantly different from the negative standard.

Sensitivity

Fish: 0.05 ppb

Shrimp: 0.05 ppb

Pork: 0.05 ppb

Chicken: 0.05 ppb

Egg: 0.15 ppb

Milk: 0.3 ppb

Specificity

The specificity of CD Florfenicol: & Thaimphenicol: ELISA kits has been determined by analyzing the cross-reactivity of the corresponding substances.

Florfenicol: 100%

Thiamphenicol: 110%

Chloramphenicol: <1%

Florfenicol amine: <1%