



User's Manual

Azaperone-Azaperol ELISA kit



DEIA-XY23



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The azaperone-azaperol ELISA is used for the quantitative analysis of azaperol and azaperone in tissue (muscle, liver, kidney) and urine samples.

Storage

2-8°C

Detection Limit

0.2 ppb for urine, 3 ppb for tissue, 5 ppb for liver and 10 ppb for kidney

Specificity

The kit uses a polyclonal rabbit antibody which shows the following cross-reactivities: azaperone (100%) and azaperol (100%).