



User's Manual

Diethylstilbestrol ELISA kit



DEIA-XY27



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The DES ELISA is used for the quantitative analysis of DES in cattle urine or tissue samples (muscle).

Storage

2-8°C

Detection Limit

0.5 ppb for milk, 0.5 ppb for egg, 0.9 ppb for liver, 0.9 ppb for tissue and 0.5 ppb for urine

Specificity

The kit uses a polyclonal rabbit antibody which shows the following cross-reactivities: DES (100%), dienestrol (10%) and hexestrol (10%).