



User's Manual

Zeranol ELISA kit



DEIA-XY47



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA kit is designed to detect zeranol in urine samples.

Storage

2-8°C

Detection Limit

0.1 ng/mL in the urine matrix

Specificity

The kit shows the following cross-reactivities: zeranol (100%), taleranol (133%), zearalanone (100%), α -zearalenol (90%), zearalenone (75%) and β -zearalenol (70%). The stilbenes (diethylstilbestrol, hexestrol and dienestrol), 17 β -oestradiol, progesterone and testosterone show a cross-reactivity < 0.01%.