



User's Manual

α -Amylase ELISA kit



DEIA-XY51



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Salivary Alpha-Amylase Immunoassay Kit (Kinetic Reaction) is specifically designed and validated for the kinetic measurement of salivary α -amylase activity. It is not intended for diagnostic use. It is intended only for research use in humans and some animals. This kit has not been validated for serum or plasma samples.

Storage

All unopened components of this kit are stable at 2-8°C until the kit's expiration date.

Precision

Intra-Assay Precision: 2.5%-7.2%

Inter-Assay Precision: 3.6%-5.8%

Sensitivity

N/A