



User's Manual

Uric Acid ELISA kit



DEIA-XY60



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Salivary Uric Acid Assay Kit was specifically designed to standardize the determination of uric acid concentrations in saliva samples. It is not intended for diagnostic use.

Storage

All unopened components of this kit are stable at 2-8°C until the kit's expiration date.

Precision

Intra-Assay Precision: 1.6%-2.5%

Inter-Assay Precision: 3.5%-6.8%

Detection Range

5 mg/dL (NIST certified)

Sensitivity

0.07 mg/dL LOD