



## User's Manual

# Acetaminophen ELISA Kit



DEIA-XYZ191



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This kit is a qualitative one-step kit designed for use as a screening device for the detection of Acetaminophen and/or their metabolites in human urine, blood or oral fluid. The kit was designed for screening purposes and is intended for forensic use only.

### Sensitivity

Compound I-50 in EIA Buffer

Acetaminophen 402 ng/mL

### Specificity

Compound Cross-Reactivity

Acetaminophen 100%

Procainamide 2.0%

Ethyl-p-amino-benzoate 1.0%

Oxyphenbutazone 0.7%

Penicillin G-Procaine 0.3%

Procaine 0.1%

Methylene Blue 0.09%