




Dextromethorphan ELISA Kit



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)

441 (Europe) Fax: 1-631-938-8221 Web: www.creative-diagnostics.com

Email: info@creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

CD's Dextromethorphan ELISA kit is a qualitative one-step kit designed for use as a screening device for the detection of trace quantities of Dextromethorphan and/or other metabolites in human urine, blood or oral fluid. The DMX drug test kit was designed for screening purposes and is intended for forensic use only.

Sensitivity

Compound I-50 in EIA Buffer

Dextromethorphan 1.1 ng/mL

Specificity

Compound Cross-Reactivity

Methoxymorphinan 118%

Dextromethorphan 100%

Hydroxymorphinan 1.7%

Dextrorphan 1.6%

Promazine 0.05%

Trimipramine 0.04%

Orphenadrine 0.04%

Imipramine 0.02%

Nortiptyline 0.02%

Amitriptyline 0.02%

Acepromazine 0.01%

Meperidine 0.01%

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

Fax: 1-631-938-8221