



User's Manual

Norepinephrine ELISA Fast Track Kit



DEIA05751



96T






This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Enzyme Immunoassay for the quantitative determination of Noradrenaline (Norepinephrine) in plasma and urine.

Storage

Store the kits at 2-8°C in a refrigerator. For more detailed information, please download the following document on our website.

Sensitivity

Analytical Sensitivity (Limit of Detection)

Urine: 1.5 ng/mL

Plasma: 50 pg/mL