



User's Manual

Phosphatidic Acid Antibody IgG/IgM ELISA Kit



DEIA1818



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Anti-Phosphatidic Acid is an indirect solid phase enzyme immunoassay (ELISA) for measurement of IgG and IgM class autoantibodies against Phosphatidic Acid in human serum or plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

IgG 0-100 GPL U/mL,

IgM: 0-100 MPL U/mL