



## User's Manual

# 1,25(OH)<sub>2</sub> Vitamin D ELISA Kit



DEIA1974V2



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Immunoenzymetric assay for the precise measurement of 1,25(OH)<sub>2</sub> Vitamin D concentrations. For research use only, not for use in diagnostic procedures.

### Storage

Before opening or reconstitution, all kits components are stable until the expiry date, indicated on the label, if kept at 2 to 8°C; except the cartridges which must be stored at room temperature.

After their reconstitution, calibrators and controls are stable for one week at 2-8°C. For longer storage periods, aliquots should be made and kept at -20°C for 2 months maximum. Avoid subsequent freeze-thaw cycles.

### Precision

Intra-assay Precision: 11.3-5.1%

Inter-assay Precision: 11.9-8.0%

### Sensitivity

1.7pg/mL.