



User's Manual

Porcine Insulin High Range ELISA Kit



DEIA2169



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The product provides a method for the direct quantitative determination of high porcine insulin levels in serum, plasma and cell cultures, without the need for high dilution of sample.

Storage

Store all contents at 2-8°C.

Precision

Inter-Lot Precision: 3.8-5.1%

Intra-assay Precision: 3.6-4.9%

Inter-assay Precision: 1.8-2.2%

Detection Range

2.5-75 ug/L

Sensitivity

0.313 µg/L