



User's Manual

17-α-OH Progesterone ELISA Kit



DEIA2233



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)

41 (Europe) Fax: 1-631-938-8221

Web: www.creative-diagnostics.com

Email: info@creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The 17- α -OH Progesterone ELISA is an enzyme immunoassay for the quantitative measurement of 17- α -OH Progesterone (17- α -OHP) in serum.

Storage

Store the kit at 2-8°C. Keep microwells sealed in a dry bag with desiccants. For more detailed information, please download the following document on our website.

Detection Range

0.034-20 ng/mL.

Sensitivity

0.034 ng/mL

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221

Email: info@creative-diagnostics.com