



User's Manual

Human Vitamin D BP ELISA Kit

REF

DEIA2702



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

For the quantitative determination of human Vitamin D Binding Protein (Vitamin D BP) concentrations in cell culture supernates, serum, plasma, saliva, urine, and human milk.

Storage

Store at 2-8°C. Do not use past kit expiration date. For more detailed information, please download the following document on our website.

Detection Range

0.15-3.74 ng/mL

Detection Limit

0.65 ng/mL