



User's Manual

Mouse Sonic Hedgehog/Shh N-Terminus ELISA Kit



DEIA2938



96T







This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

For the quantitative determination of mouse Sonic Hedgehog N-Terminus (Shh-N) concentrations in cell culture supernates, mouse serum, and plasma.

Storage

Store at 4°C for frequent use, at -20°C for infrequent use. Avoid multiple freeze-thaw cycles.

Return unused wells to the foil pouch containing the desiccant pack, reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.

Precision

Intra-assay Precision:
2.4-5.0%

Inter-assay Precision:
4.9-10.4%

Detection Range

0.032-2.37 pg/mL

Sensitivity

1.22 pg/mL