



User's Manual

Measles IgA ELISA Kit



DEIA360



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Measles IgA Antibody ELISA Test Kit has been designed for the the detection and the quantitative determination of specific IgA antibodies against Measles in serum and plasma.

Storage

For more detailed information, please download the following document on our website.

Precision

Intra-assay-Precision: 12.5%

Inter-assay-Precision: 6.9-15.8%

Inter-Lot-Precision: 2.7-12.2%

Sensitivity

1.35 U/mL