



User's Manual

# Human Prorenin ELISA Kit



DEIA4768



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Human Prorenin assay is intended for the quantitative determination of Prorenin in biological fluids. Active renin will not be detected by this assay. Prorenin is measured directly by ELISA without pretreatment of samples or conversion to renin.

### Storage

All kit components must be stored at 4°C. For more detailed information, please download the following document on our website.

### Detection Range

The assay measures human Prorenin in the 0.01-10 ng/mL range.