



User's Manual

Prostaglandin F2 α ELISA Kit



DEIA5002



96T, 480T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay is based on the competition between PGF2 α and a PGF2 α -acetylcholinesterase (AChE) conjugate (PGF2 α Tracer) for a limited amount of PGF2 α -specific rabbit antiserum.

Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.

Detection Range

3.9-500 ng/mL

Detection Limit

3.9 ng/mL