




User's Manual

Enterolactone EIA Kit



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221

PRODUCT INFORMATION

Intended Use

CD's Enterolactone EIA Kit is a competitive assay that can be used for quantification of enterolactone in plasma, urine, and other sample matrices.

Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.

Detection Limit

70 pg/mL

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221

Email: info@creative-diagnostics.com