




User's Manual

11-dehydro Thromboxane B2 EIA Kit



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe) Fax: 1-631-938-8221

PRODUCT INFORMATION

Intended Use

CD's 11-dehydro TXB2 EIA Kit is a competitive assay that can be used for quantification of 11-dehydro TXB2 in urine and other sample matrices. The EIA typically displays an IC50 (50% B/B0) value of 120 pg/mL and a detection limit (80% B/B0) of approximately 34 pg/mL.

Storage

This kit will perform as specified if stored as directed at -20°C and used before the expiration date indicated on the outside of the box. For more detailed information, please download the following document on our website.

Detection Limit

34 pg/mL

Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)

Fax: 1-631-938-8221

Email: info@creative-diagnostics.com