



User's Manual

Total sAPP high sensitive ELISA Kit



DEIA6193



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay kit is capable for the quantitative determination human total sAPP in serum, EDTA plasma, cerebrospinal fluid and cell culture supernatant as a total amount.

Storage

Storage Condition: 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0.39-25 ng/mL

Sensitivity

0.06 ng/mL

The sensitivity for this kit was determined using the guidelines under the National Committee for Clinical Laboratory Standards (NCCLS) Evaluation Protocols. (National Committee for Clinical Laboratory Standards Evaluation Protocols, SC1, (1989) Villanova, PA: NCCLS.)