



User's Manual

Androstane-diol Glucuronide ELISA Kit

REF

DEIA6239



96T



RUO


This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Androstane-diol Glucuronide ELISA is for the direct quantitative determination of 5 α -Androstane-3 α , 17 β -diol Glucuronide (3 α Diol G) by enzyme immunoassay in human serum.

Storage

The kit is shipped at ambient temperature and should be stored at 2-8°C. Keep away from heat or direct sun light. The storage and stability of specimen and prepared reagents is stated in the corresponding chapters. Once opened, Standards and Controls should be used within 14 days or aliquoted and stored frozen. Avoid multiple freezing and thawing cycles.

The Microtiter Plate and other reagents are stable up to 12 months or indicated on label, when stored at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0-50 ng/mL

Sensitivity

0.1 ng/mL