



## User's Manual

# Bovine MALACT ELISA Kit



DEIA8721



10 x 100 tests



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Bovine MALACT ELISA Kit is for the detection of Bovine Alpha-Lactalbumin in milk. The Set contains sufficient components to perform 1000 single well assays.

The set performance has been optimized for the stated protocol using the materials listed and standard dilutions from 50-0.78 ng/mL of Bovine Alpha-Lactalbumin. For alternative assay conditions, the operator must determine appropriate dilutions of reagents.

### Storage

Store all reagents at 2-8°C

### Detection Range

Assay Range:

0.78-50 ng/mL

Suggested standard curve points are 50 ng/mL, 25 ng/mL, 12.5 ng/mL, 6.25 ng/mL, 3.125 ng/mL, 1.56 ng/mL, 0.78 ng/mL and 0 ng/mL (blank).