



User's Manual

Bovine LTF ELISA Kit



DEIA8739



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Bovine LTF ELISA Kit is intended for the detection of Bovine Lactoferrin in milk, colostrum, or plasma. Other biological fluids that contain Bovine Lactoferrin, such as urine, tears, or saliva, may be suitable samples.

Storage

Store all reagents at 2-8°C

Detection Range

Assay Range:

7.8-500 ng/mL

Suggested standard curve points are 500 ng/mL, 250 ng/mL, 125 ng/mL, 62.5 ng/mL, 31.25 ng/mL, 15.625 ng/mL, 7.8 ng/mL and 0 ng/mL (blank).