



## User's Manual

# Prostaglandin E2 ELISA Kit



DEIABL275



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This competitive ELISA kit is for determination of PGE2 (Prostaglandin E2) levels in biological samples.

Each kit can be used for triplicate analyses of up to 24 samples contains using a 96 well plate format, and contains a vial of PGE2 standard, a vial of PGE2-conjugated horseradish peroxidase (HRP), and buffers for sample and HRP dilutions, and plate washing.

### Storage

This kit will obtain optimal results if all of the components are stored at the proper temperature prior to use. Items should be stored at the designated temperatures upon receipt of this kit.

All components are stored below -20°C and should not be re-frozen and thawed more than necessary.