



Fluorescent Brightener Ksn [KLH] (DAGA-109K)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature S	Synthetic
Expression System N	J/A
Conjugate K	(LH
Applications In	mmunogen
Procedure N	None
Format L	iquid
Concentration B	Batch dependent - please inquire should you have specific requirements.
Size 1	mg
Buffer 0	0.01M pH7.4 PBS
Preservative N	None
Storage S	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Keywords Fluorescent Brightener KS-N;Fluorescent Brightener KSN

45-1 Ramsey Road, Shirley, NY 11967, USA

Email:info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221