



Native Campylobacter Coli (DAGA-3022)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Nature	Native
Expression System	Campylobacter coli clinical strain obtained
Conjugate	Unconjugated
Applications	Suitable for use as control in Lateral Flow assay and as an immunogen. May be used to spike stool samples. Each laboratory should determine an optimum working titer for use in its particular application. Other applications have not been tested but use in such assays should not necessarily be excluded.
Format	Lysate
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	0.1 ml
Buffer	137 mM Sodium Chloride (NaCl), 2.7 mM Potassium Chloride (KCl), 10 mM Disodium Phosphate (Na ₂ HPO ₄), 1.8 mM Dihydrogen Phosphate (H ₂ PO ₄), pH 7.4 + 0.2.
Preservative	0.09% Sodium Azide
Storage	Store at 2-8°C for < 1 week. Long term storage at -20°C. Aliquot to avoid multiple freeze/thaw cycles.