



Native Salmonella typhimurium antigen (DAGA-3051)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	ELISA, Lateral Flow
Format	Liquid; Inactivated Native Antigen
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	100 μL
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservative	0.09% Sodium Azide
Storage	Short Term (< 1 week): 2-8 centigrade. Long Term: -20 centigrade. Avoid repeated freezing and thawing.

Email: info@creative-diagnostics.com

Tel: 1-631-624-4882 Fax: 1-631-938-8221