



Salmonella pullorum and gallinarum positive serum (DAGF-121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

or red blood cells) nor a clotting factor; it is the blood plasma not including the fibrinogens. Serum includes all proteins not used in blood clotting (coagulation) and all the electrolytes,

antibodies, antigens, hormones, and any exogenous substances (e.g., drugs and

microorganisms). A study of serum is serology, and may also include proteomics. Serum is

used in numerous diagnostic tests, as well as blood typing.

Conjugate Unco	njugated
-----------------------	----------

Format Liquid

Concentration Batch dependent - please inquire should you have specific requirements.

Size 1 ml

Preservative None

Storage Store at -20°C or lower. Avoid repeated freeze/thaw cycles

BACKGROUND

Keywords S. pullorum; S. gallinarum; salmonella pullorum; salmonella gallinarum; serum