



Bacterial Vaginosis Rapid Test Kit (Sialidase) (DST-CL003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Intended Use	Bacterial Vaginosis Rapid Test Kit was designed for detection of abnormal sialidase activity of vaginal discharge of women suffering from BV.
Storage	The test kit can be stored at temperature (4°C to 26°C) in the sealed pouch for 2 years. The test kits should be kept away from direct sunlight, moisture and heat.
Precision	The intra and inter-batch imprecision of reference have no difference.
Sensitivity	7.8μ/L