



Clenbuterol Residue Rapid Test Strip (Serum) (DTS002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Size	40T
Intended Use	Clenbuterol device is for rapid test to qualitatively detect the Clenbuterol in Serum sample at the sensitivity of 3 µg/kg. It only takes approx. 5–8 min.
Storage	Store at 15-25°C, DO NOT FREEZE or use beyond the expiration date. The shelf life is 12 months.
Sensitivity	To acquire the exact sensitivity, reduplicative experiment has been done on the sample containing 3 µg/kg Clenbuterol.